


Zinc oxide adhesive plaster



- Ref no:1.11.1
- Material:cotton
- natural rubber and zinc oxide adhesive
- easy to take
- Size: 5cm*4.5m,7.5cm*4.5m,10cm*4.5m
- Packing can be customized

			
Zinkoxid Band	Bande d'oxyde de zinc	Cinta de óxido de zinc	酸化亜鉛テープ

Cotton elastic bandage



- Ref no:1.11.2
- Material: 85%cotton+spandex
- natural rubber and zinc oxide adhesive
- easy to take
- Size: 5cm*4.5m,7.5cm*4.5m,10cm*4.5m
- Packing can be customized

			
Baumwolle elastische Binde	85%cotton+spandex	Vendaje elástico de algodón	コットン弾性包帯

High Elastic Bandage



- Ref no:1.11.3
- Perfect for sports medicine and sports-related injuries
- Value and durability make it ideal for everyday-use
- Strong enough to withstand laundering
- Size: 5cm*4.5m,7.5cm*4.5m,10cm*4.5m
- Packing can be customized

			
hohe elastische Binde	haute bandage élastique	alta vendaje elástico	高弾性包帯

Cotton Crepe Bandage



- Ref no:1.11.4
- cotton material
- meet B.P standard
- With or without blue/red lines
- Moderate elasticity
- Hospital quality
- Size: 5cm*4.5m,7.5cm*4.5m,10cm*4.5m
- Packing can be customized

			
Baumwolle elastische Binde	crépon de coton bandage	crepé de algodón vendaje	綿クレープ包帯

Heavy Elastic Bandage



- Ref no:1.11.5
- Heavy and stable stretch
- Constant unwinding tension
- Strong and reliable stickness
- Size: 5cm*4.5m,7.5cm*4.5m,10cm*4.5m
- Packing can be customized

			
PE Schürze	PE Tablier	Pesado vendaje elástico	ヘビー弾性包帯

Self Closure Bandage



- Ref no:1.11.6
- Velcro both ends.Requires no pins,clips or tape
- 52%cotton and 20% spandex and 38% latex
- Available in multiple package, Normal packing for individual is flow wrapped
- Size: 5cm*4.5m,7.5cm*4.5m,10cm*4.5m
- Packing can be customized

			
Selbst Closure Bandage	Self Closure Bandage	Cierre auto vendaje	自己閉鎖包帯

Self adhesive Elastic Bandage



- Ref no:1.11.7.1
- Material: cotton or non woven
- Non-sticky to hair or skin, Hand tear
- Provide controlled compression
- Soft, water resistant
- Size: 5cm*4.5m, 7.5cm*4.5m, 10cm*4.5m
- Packing can be customized

			
PE Schürze	Auto adhésif Bandage élastique	auto-adhesivo vendaje elástico	自己接着剤弾性包帯

Stockinet bandage(100%Polyester)



- Ref no:1.11.7.2
- Material: 100%Polyester
- Provide controlled compression
- Soft, water resistant
- Size: 5cm*20m, 7.5cm*20m, 10cm*20m
- Packing can be customized

			
stockinet spielen Big Brother	stockinet jouer Big Brother	Telas de jugar a Gran Hermano	メリヤスビッグブラザーを 再生

Stockinet bandage(100%Cotton)



- Ref no:1.11.7.3
- Material: 100%Cotton
- Provide controlled compression
- Soft, water resistant
- Size: 5cm*20m, 7.5cm*20m, 10cm*20m
- Packing can be customized

			
stockinet spielen Big Brother	stockinet jouer Big Brother	Telas de jugar a Gran Hermano	メリヤスビッグブラザーを 再生

Conforming Pressure Bandages



- Ref no:1.11.8
- conforming bandage and non-stick pad
- Soft and comfortable
- No stimulation to the skin, not allergic
- Size: 6cm*4m, 8cm*4m, 10cm*4m
- Packing can be customized

			
Konformer Druck Bandagen	Conformes Bandages de pression	Los vendajes de presión conforme	圧力包帯を従わせます

Tubular elastic bandage



- Ref no:1.11.9.1
- Breathable and comfortable to be worn in shoes and underclothing
- Washable
- 5cm*4.5m, 7.5cm*4.5m, 10cm*4.5m
- Packing can be customized

			
Tubular elastische Binde	bandage élastique tubulaire	venda elástica tubular	管状弾性包帯

Tubular bandage



- Ref no:1.11.9.2
- Breathable and comfortable to be worn in shoes and underclothing
- Washable
- 6.25*20m, 7.5cm*20m, 8.75cm*20m
- Packing can be customized

			
Tubular Binde	bandage tubulaire	venda tubular	管状包帯

Tubular Net Bandages



- Ref no:1.11.10
- Easy wound inspection and redressing.
- Washable and reusable.
- Size: 5cm*4.5m, 7.5cm*4.5m, 10cm*4.5m
- Packing can be customized

			
Netzschlauchverband	Tubular Net Bandage	Vendaje tubular neto	筒状ネット包帯

Plaster of Paris Padding



- Ref no:1.11.11
- Made by Smooth POP Powder,
- Low Powder loss.
- Short Immersion time.
- Smoother and creamy finish, gives early, easy and excellent mouldability and higher stability.
- Size:5cm*270cm,7.5cm*270cm,10cm*270cm
- Packing can be customized

			
Plaster of Paris Padding	Le plâtre de Paris Rembourrage	Yeso de París Relleno	パリのパディングの石膏

Gauze bandage



- Ref no:1.11.12
- Material:100%bleached cotton
- Yarn: 40S 24*20
- Good quality, high liquid absorbency.
- Short Immersion time.
- For disposable use
- Size:5cm*9.1m,7.5cm*9.1m
- Packing can be customized

			
Mullbinden	bandages de gaze	vendas de gasa	ガーゼの包帯

Undercast padding



- Ref no:1.11.13
- Material:cotton+viscose ,polyester+viscose
- Size:7.5cm*2.7m,10cm*2.7cm or available
- Weight:90gsm
- Highly conformable
- Easily tearable
- High soft
- Packing can be customized

			
under Polsterung	padding undercast	acolchado undercast	アンダーパディング