


## Pillow case



- Ref no. : 1.1.6.1.1
- Material: PP Nonwoven
- Edge unfolded
- Size: 50\*70cm size can be customized
- Packing: 10pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                   |                                                                                   |                                                                                   |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| Kissenbezug                                                                       | Taie d'oreiller                                                                   | funda de almohada                                                                 | ピローケース                                                                              |

## Pillow case



- Ref no. : 1.1.6.1.2
- Material: PP Nonwoven
- Edge folded
- Size: 50\*70cm size can be customized
- Packing: 10pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Kissenbezug                                                                         | Taie d'oreiller                                                                     | funda de almohada                                                                   | ピローケース                                                                                |

## PP bed sheet



- Ref no.:1.1.6.2.1
- Material: PP Nonwoven
- 100% new material
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| PP Bettlaken                                                                        | Drap PP                                                                             | Sábana PP                                                                           | PP のベッドシーツ                                                                            |

## PP bed sheet



- Ref no.:1.1.6.2.2
- Material: PP Nonwoven
- Normal material with recycled material
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                   |                                                                                   |                                                                                    |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| PP Bettlaken                                                                      | Drap PP                                                                           | Sábana PP                                                                          | PP のベッドシート                                                                          |

## SMS bed sheet



- Ref no.:1.1.6.2.3
- Material: SMS material
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| SMS Bettlaken                                                                       | SMS de base d t Shee                                                                | Sábana SMS                                                                           | SMS ベッドシート                                                                            |

## PP+PE bed sheet



- Ref no.:1.1.6.2.4
- Material: PP+PE material
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Tsu p + Bae Basis d Shee t                                                          | Tsu p + base Bae d Shee t                                                           | Sábana PP + PE                                                                       | PP+ PE ベッドシート                                                                         |

## Impregnated bed sheet



- Ref no.:1.1.6.2.5
- Material: Impregnated cloth
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                   |                                                                                   |                                                                                   |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| Es mp Re g a in d Basis d<br>Shee t                                               | Imprégnés Bed Sheet                                                               | Hay pf Re g en una base de d<br>d t Shee                                          | 含浸されたベッドシート                                                                         |

## Bed sheet with antistatic treatment



- Ref no.:1.1.6.2.6
- Material: PP Nonwoven
- With antistatic
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Bettlaken mit antistatischem                                                        | Drap avec antistatique                                                              | Sábana con antiestático                                                             | 帯電防止とベッドシート                                                                           |

## Bed sheet with hydrophilic treatment



- Ref no.:1.1.6.2.7
- Material: PP Nonwoven
- With hydrophilic
- Size: 80\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Bettlaken mit hydrophilen                                                           | Base de d Shee t Wie e                                                              | Base d Shee t Wie                                                                   |                                                                                       |

## 3-ply bed sheet



- Ref no.:1.1.6.2.8
- Material: 3-ply PP material
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                   |                                                                                   |                                                                                   |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| 3-lagig Bettlaken                                                                 | Drap de lit 3 plis                                                                | Sábana de 3 capas                                                                 | 3 プライベッドシート                                                                         |

## 4-ply bed sheet



- Ref no.:1.1.6.2.9
- Material: 4-ply PP material
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| 4-lagig Bettlaken                                                                   | Drap de lit 4 plis                                                                  | 4 sobre capa base d t Shee                                                          | 4 プライのベッドシート                                                                          |

## PP bed cover



- Ref no.:1.1.6.3.1
- Material: PP non woven material
- With elastic in 4 corner, with elastic in both sides or elastic around sides
- Size: 70\*180cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Tsu p-Basis d Kove r                                                                | Couvre-lit PP                                                                       | Tsu base de p r d Kovell                                                            | PP のベッドカバー                                                                            |

## PP+PE bed cover



- Ref no.:1.1.6.3.2
- Material: PP+PE material
- With elastic in 4 corner or with elastic in both sides or elastic around sides
- Size: 70\*180cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                   |                                                                                   |                                                                                   |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| PP + PE Bettdecke                                                                 | Couverture de lit PP + PE                                                         | Cubierta espacio de PP + PE                                                       | PP+ PE のベッドカバー                                                                      |

## Non-woven blanket



- Ref no.:1.1.6.4.1
- Material: PP nonwoven and with filler of cotton
- Size: 70\*200cm size can be customized
- Packing: 10pc/PE bag, 5bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Die Non Ven b Anke t                                                                | Le non Ven b Anke t                                                                 | No tejido manta                                                                     | 不織布毛布                                                                                 |

## Pillow



- Ref no.:1.1.6.5.1
- Material: PP Nonwoven
- All sides heated by ultrasonic or sewed by machine
- Size: 50\*50cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized

|                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Kissen                                                                              | oreiller                                                                            | Pi'o w                                                                               | 枕                                                                                     |

## PP+RECYCLE PAPER Blanket



- Ref no.:1.1.6.6.1

2-ply PP20g/sm with 4-ply recycle paper 20g/sm

Size: 110\*190cm

Colou: white

Total: 250g/pc, 6-layers/pc

Packing: 1pc/bag, 75bags/ctn, 60\*40\*40cm

Colour and size could be changed by customer's deamnd.

Anti-static and Flame resistance test passed

|                                                                                   |                                                                                   |                                                                                    |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| Decke mit PP + RECYCLE<br>PAPER                                                   | Couverture avec PP +<br>RECYCLE PAPIER                                            | Manta con papel PP +<br>RECICLAR                                                   | PP の+の再生紙とブランケ<br>ット                                                                |

## PP+ Needle-punched nonwoven fabric Blanket



Ref no.:1.1.6.6.2

-Material:2-ply PP with Needle-punched nonwoven fabric

-Size:110\*190cm,110\*130,130\*210,150\*190-Purpose:keep warm, for refugee

-Packing: 1pc/bag,20 or 40bags/ctn

- Process:sewing/ ultrasonic

-Color option: white,green,blue

Anti-static and Flame resistance test passed

|                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| PP Nadelvlies Stoff Blanket                                                         | PP + aiguilleté non tissé<br>Couverture en tissu                                    | PP + punzonado manta de<br>tela no tejida.                                           | PP ニードルパンチ不織布ブ<br>ランケット                                                               |

## PP+Polyester Blanket

Ref no.:1.1.6.6.3



- Material: 2-ply PP with Polyester fabric
- Size:110\*190cm,110\*130,130\*210,150\*190
- Purpose:keep warm, for refugee
- Packing: 1pc/bag,20 or 40bags/ctn
- Process:sewing/ ultrasonic
- Color option: white,green,blue
- Anti-static and Flame resistance test passed

|                                                                                   |                                                                                   |                                                                                    |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| PP + Polyester-Decke                                                              | Tsu p + Poe s sur r B Anke t                                                      | Tsu p + Poe s de r B Anke t                                                        | PP+ポリエステル毛布                                                                         |

## Protective sheet



- Ref no: 1.1.6.7
- It is made of 1ply of paper and 1 ply of PE film.
- Specification: 80\*140cm, 80\*170cm, 80\*210cm or customized
- Color: white or blue
- Packed in sheet or in roll
- Packing:25pcs/bag,4bags/ctn

|                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Schutzfolie                                                                         | feuille de protection                                                               | lámina protectora                                                                    | 保護シート                                                                                 |

## Field drape



- Ref no: 1.1.6.8
- It is made from paper and PE film.
- 1-ply of paper+1ply pe film pr 2ply of paper +1ply pe film
- And with strengthening stripes in between.
- Specification: 80\*140cm, 80\*170cm, 80\*210cm
- The strengthening strips can be 8pcs, 24pcs or 48pcs.
- Packing:25pcs/bag,4bags/ctn

|                                                                                     |                                                                                     |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| Schutzfolie                                                                         | feuille de protection                                                               | lámina protectora                                                                   | 保護シート                                                                                 |

## Exam paper roll



- Ref no: 1.1.6.9
- It is made of wet strong paper
- 1ply of paper or two ply of paper
- Specification: 18"x225' , 21"x225' or customized
- Packing:12rolls/ctn

|                                                                                   |                                                                                   |                                                                                    |                                                                                     |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| Schutzfolie                                                                       | feuille de protection                                                             | lámina protectora                                                                  | 保護シート                                                                               |

## CPE bed cover



- Ref no:1.1.6.10
- It is made of CPE material
- widely applied in hospital, clinics, spa room and massage room
- Keep cleaning
- Size: 20\*90\*200cm
- Packing:100pcs/ctn

|                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |
| CPE Bettdecken                                                                      | couvre-lits CPE                                                                     | colchas CPE                                                                          | CPE のベッドカバー                                                                           |