

Isolation gown with elastic cuff



- Ref no. : 1.1.3.4.1
- Material: PP Nonwoven or SMS and so on
- With elastic cuff neck and 2 waist tape
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

			
Isolation Kleid	robe d'isolement	bata de aislamiento	アイソレーションガウン

Isolation gown with knitted cuff



- Ref no. : 1.1.3.4.2
- Material: PP Nonwoven or SMS and so on
- With knitted cuff neck and 2 waist tape
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

			
Isolation Kleid	robe d'isolement	bata de aislamiento	アイソレーションガウン

Surgical gown



- Ref no. : 1.1.3.7.1
- Material: SMS
- With knitted cuff Velcro neck and 2 waist tape
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

			
OP-Kittel	blouse chirurgicale	bata quirúrgica	手術衣

Reinforced Isolation gown



- Ref no. : 1.1.3.4.3
- Material: PP+PE
- With knitted cuff, Velcro neck and waist tape.
- Reinforced with PP+PE for the whole gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized



OP-Kittel verstärkt



blouse chirurgicale renforcée



bata quirúrgica reforzada



手術用ガウンを強化します

Reinforced Isolation gown



-Ref no. : 1.1.3.4.4

- Material: PP+PE
- With elastic cuff, Velcro neck and waist tape.
- Reinforced with PP+PE for sleeve and front area.
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized



OP-Kittel verstärkt



blouse chirurgicale renforcée



bata quirúrgica reforzada



手術用ガウンを強化します

Reinforced Surgical gown



- Ref no. : 1.1.3.7.3
- Material: SMS
- With knitted cuff, Velcro neck and 4 waist tape.
- Reinforced with PP+PE in front of gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized



OP-Kittel verstärkt



blouse chirurgicale renforcée



bata quirúrgica reforzada



手術用ガウンを強化します

Reinforced Isolation gown by ultrasonic



- Ref no. : 1.1.3.4.3.1
- Material: PP+PE
- With elastic cuff, neck and waist tape. by ultrasonic.
- Reinforced with PP+PE for the whole gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized



OP-Kittel verstärkt



blouse chirurgicale renforcée



bata quirúrgica reforzada



手術用ガウンを強化します

Surgical gown by ultrasonic



- Ref no. : 1.1.3.7.1.1
- Material: SMS
- With knitted cuff Velcro neck and 2 waist tape, by ultrasonic.
- Size: 110*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized



OP-Kittel verstärkt



blouse chirurgicale renforcée



bata quirúrgica reforzada



手術用ガウンを強化します

Reinforced Surgical gown by ultrasonic



- Ref no. : 1.1.3.7.3.1
- Material: SMS
- With knitted cuff, Velcro neck and 4 waist tape, by ultrasonic.
- Reinforced with PP+PE in front of gown
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white, blue or all color can be customized



OP-Kittel verstärkt



blouse chirurgicale renforcée



bata quirúrgica reforzada



手術用ガウンを強化します

Patient gown



- Ref no. : 1.1.3.4.1.2
- Material: PP Nonwoven or SMS and so on
- without sleeve neck and 2 waist tape
- Size: 120*140cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

			
Patientenhemd	robe du patient	bata de paciente	患者ガウン

Childbirth shirt



- Ref no. : 1.1.3.10
- Material: PP Nonwoven
- Size: 180*110cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

			
Chi Id Schönheit rth auf rt	La beauté de rth de Chi à la température ambiante	Chi Id rth belleza a rt	出産のシャツ

Patient exam capes



- Ref no. : 1.1.3.10.1
- Material: PP Nonwoven
- No sleeve, neck and waist with tie
- Size: 75*52cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: white, blue or all color can be customized

			
Prüfung Umhänge	capes d'examen	cabos del examen	試験の岬

Spunlace surgical gown



- Ref no.:1.1.3.7.4
- Material: Spun lace wood pulp
- With knitted cuff, Velcro neck and 4 waist tape
- Size: 110*120cm, size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: blue or all color can be customized

			
Spunlaced-OP-Kittel	dentelle filée blouse chirurgicale	bata quirúrgica de encaje hilado	スパンレース手術用ガウン

Scrub suits



- Ref no.:1.1.3.8
- Material: PP Nonwoven or SMS
- Short sleeve With V neck or round neck With or without pocket, Elastic or tie waist
- Top Size: 70*110cm pant size: 105*100cm, size can be customized
- Packing: 1pcs/PE bag, 50bags/ctn, can be customized
- Color option: blue or all color can be customized

			
Bereichskleidung	costumes gommage	trajes de matorral	スクラブスーツ

Exam Gown



- Ref no.:1.1.3.9
- Material: PP Nonwoven or SMS
- Short sleeve or no sleeve with or without pocket,
- Top Size: 110*120cm, size can be customized
- Packing: 1pcs/PE bag, 50bags/ctn, can be customized
- Color option: blue or all color can be customized

			
Prüfung Kleid	robe examen	bata de examen	試験のガウン

Coverall



- Ref no.:1.1.3.3
- Material: PP Non woven or SMS
- With hood or without hood, zipper, elastic cuff and bottom, waist
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

			
Overall	coverall	sobretudo	カバーオール

Coverall Type 5/6



- Ref no.:1.1.3.3.3
- Material: PP Non woven + Breathable film
- With hood, zipper, elastic cuff and bottom, waist
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or blue, color can be customized

			
Overall	coverall	sobretudo	カバーオール

Coverall Type4/5/6



- Ref no.:1.1.3.3.4
- Material: PP Non woven + Breathable film
- Elastic in waist area providing body-fitting of a coverall, elastic in cuffs, legs and hood, loose shape with gore in groin area, three-part hood , self-adhesive cover of a zipper, adhesive chin cover, with elastic thumb loops
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or blue, color can be customize

			
Overall	coverall	sobretudo	カバーオール

Lab coat



-Ref no.:1.1.3.5

Material: PP Non woven or SMS

- With collar, Elastic or knitted cuff Velcro or snap closure
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

			
Laborkittel	blouse de laboratoire	bata de laboratorio	白衣

Lab coat



-Ref no.:1.1.3.5.1

Material: PP Non woven or SMS

- With Knitted collar, Elastic or knitted cuff Velcro or snap closure
- Size: 120*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

			
Laborkittel	blouse de laboratoire	bata de laboratorio	白衣

Non woven Kimono



-Ref no.:1.1.3.14

-Material: PP Non woven

- With one tie With one pocket in short sleeve
- Size: 100*140cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

			
non-woven Kimono	kimono non tissé	kimono no tejidas	不織布着物

Non woven shirt



- Ref no.:1.1.3.11
- Material: PP Non woven/ PP+PE
- short sleeve and v neck
- Size: 50*120cm size can be customized
- Packing: 10pcs/PE bag, 10bags/ctn, can be customized
- Color option: White, color can be customized

			
Die do Oh p Theorie Ven	Ven le do de la théorie de p Oh	Ven el do de Teoría p Oh	不織布エプロン

Non woven pant



- Ref no.:1.1.3.12
- Material: PP Non woven
- With elastic or tie waist Elastic bottom
- Size: 105*110cm size can be customized
- Packing: 1pc/PE bag, 50bags/ctn, can be customized
- Color option: white or black color can be customized

			
non-woven Verräter	non tissé traître	no tejida traidor	不織布裏切り者

Non woven apron



- Ref no.:1.1.3.16
- Material: PP Non woven/ PP+PE
- With tie in waist
- Size: 75*100cm size can be customized
- Packing: 50pcs/box, 8boxes/ctn, can be customized
- Color option: White, color can be customized

			
Die do Oh p Theorie Ven	Ven le do de la théorie de p Oh	Ven el do de Teoría p Oh	不織布エプロン

Non woven underwear



- Ref no.:1.1.3.1.1
- Material: PP Non woven
- With elastic
- Size: 52*25cm size can be customized
- Packing: 1pc/PE bag, 1000bags/ctn, can be customized
- Color option: white or black color can be customized

			
r r Hardware in Ja Ven die Do von	r r Hardware dans Ouais Ven le do de	La ropa interior no tejida	不織布下着

Non woven underwear



- Ref no.:1.1.3.1.2
- Material: PP Non woven
- With elastic
- Size: 44*16cm size can be customized
- Packing: 1pc/PE bag, 1000bags/ctn, can be customized
- Color option: white or black color can be customized

			
r r Hardware in Ja Ven die Do von	r r Hardware dans Ouais Ven le do de	La ropa interior no tejida	不織布下着

Non woven bra



- Ref no.:1.1.3.15
- Material: PP Non woven
- With elastic
- Size: 135*20cm, size can be customized
- Packing: 50pc/PE bag, 500bags/ctn, can be customized
- Color option: white or black color can be customized

			
Non-Woven-BH	Non tissé soutien-gorge	Sujetador no tejida	不織布ブラジャー

PE sleeves



Ref no: 1.1.3.17

Made of PE. With comfortable elastic opening

- waterproof and dustproof
- light weight, convenience, flexible, economical, soft
- Available in different thickness and color
- Size and packing can be customized

			
PE-Mantel	PE manchon	manga de PE	PE スリーブ

PE Apron



Ref no: 1.1.3.18

Made of LDPE, HDPE

- Size: 80*120cm, 80*125cm
- Color: white
- Packing: 1000pcs/ctn

			
PE Schürze	PE Tablier	PE delantal	PE エプロン

PE visit gown



-Ref no: 1.1.3.19

-Made of PE.

-Size: 120*145cm

-Color: white

-Packing: 100pcs/ctn

			
PE Schürze	PE Tablier	PE delantal	PE エプロン

CPE gown



- Ref no: 1.1.3.20
- Made of CPE.
- Size:120*198cm or other size
- Color: blue
- Packing:100pcs/ctn

			
CPE Schürze	CPE Tablier	CPE delantal	CPE エプロン

Dental bib



- Ref no: 1.1.3.21
- Size: 33*45cm
- wood pulp paper + water-proof PE backing
- 2 plies or 3plies available
- Packing: 125pcs/ bag 4bags/ctn
- Color option: different color

			
dental bib	babero dental	Den era l belleza b	歯科よだれかけ

Dental Apron Finland type



- Ref no: 1.1.3.22.1
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- with pocket in front
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental Schürze Finnland Art	Type Tablier dentaire Finlande	Tipo Delantal Dental Finlandia	歯科エプロンフィンランド式

Dental Apron Finland type



- Ref no: 1.1.3.22.2
- Made of wet strong paper and water proof PE backing
- With adhesive sticker to guarantee comfort wear
- with pocket in front
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental Schürze Finnland Art	Type Tablier dentaire Finlande	Tipo Delantal Dental Finlandia	歯科エプロンフィンランド 式


Dental Apron Italy type



- Ref no: 1.1.3.23.1
- Made of wet strong paper and water proof PE backing - With adhesive sticker to guarantee comfort wear
- With overhead
- Packing: 300pcs/ctn
- Color option: blue white pink

			
Dental Schürze Italien Art	Type Tablier dentaire Italie	Tipo Delantal Dental Italia	歯科エプロンイタリアタイプ

Dental Apron Italy type



- Ref no: 1.1.3.23.2
- Made of wet strong paper and water proof PE backing
- With adhesive sticker to guarantee comfort wear
- With pocket backside
- Packing: 300pcs/ctn
- Color option: blue white pink

			
Dental Schürze Italien Art	Type Tablier dentaire Italie	Tipo Delantal Dental Italia	歯科エプロンイタリアタイプ

Dental Apron German type



- Ref no: 1.1.3.24.1
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- With pocket backside
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental Schürze deutschen Typ	Tablier dentaire de type allemand	Delantal Dental tipo alemán	歯科エプロンドイツタイプ

Dental Apron German type



- Ref no: 1.1.3.24.2
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- With pocket backside
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental Schürze deutschen Typ	Tablier dentaire de type allemand	Delantal Dental tipo alemán	歯科エプロンドイツタイプ

Dental Apron American type



- Ref no: 1.1.3.25.1
- Made of wet strong paper and water proof PE backing
- With a smart strip to guarantee comfort wear
- With tie, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental chürze amerikanische Art	Type dentaire Tablier américain	Tipo Delantal Dental Americana	歯科エプロンアメリカンタイ プ

Dental Apron American type



- Ref no: 1.1.3.25.2
- Made of wet strong paper and water proof PE backing
- overhead, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental chürze amerikanische Art	Type dentaire Tablier américain	Tipo Delantal Dental Americana	歯科エプロンアメリカンタイ プ

Dental Apron Israel type



- Ref no: 1.1.3.26.1
- Made of wet strong paper and water proof PE backing
- overhead, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental Schürze Israel Art	Type Tablier dentaire Israël	Tipo Delantal Dental Israel	歯科エプロンイスラエルタイ プ

Dental Apron Israel type



- Ref no: 1.1.3.26.2
- Made of wet strong paper and water proof PE backing
- overhead, no pocket
- Packing: 300pcs/ctn
- Color option: blue white

			
Dental Schürze Israel Art	Type Tablier dentaire Israël	Tipo Delantal Dental Israel	歯科エプロンイスラエルタイ プ

Dental Apron Spain type



- Ref no: 1.1.3.27.1
- Made of wet strong paper and water proof PE backing
- There are two tapes on the shoulder. It's convenient.-
- Packing: 300pcs/ctn

			
Spanien Lätzchen	Espagne Bibs	España delantal dental	歯科スペインビブ

Special Dental Apron



- Ref no: 1.1.3.28.1
- Made of wet strong paper and water proof PE backing
- Size: 30*40cm overhead, no pocket
- Packing: 300pcs/ctn
- Color option: white

			
dental Schürze	tablier dentaire	delantal dental	歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.2
- Made of wet strong paper and water proof PE backing
- Size: 40*76cm,
- Packing: 300pcs/ctn
- Color option: white, blue

			
dental Schürze	tablier dentaire	delantal dental	歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.3
- Made of wet strong paper and water proof PE backing
- Size: 40*124cm,
Packing: 300pcs/ctn
- Color option: Blue, white



dental Schürze



tablier dentaire



delantal dental



歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.4
- Made of wet strong paper and water proof PE backing
- Size: 37*63cm. pocket with paper outside
Packing: 300pcs/ctn
- Color option: Blue, white



dental Schürze



tablier dentaire



delantal dental



歯科エプロン

Special Dental Apron



- Ref no: 1.1.3.28.5
- Made of wet strong paper and water proof PE backing
- Size: 37*46, 37*62cm. pocket with paper outside
Packing: 300pcs/ctn
- Color option: Blue, white



dental Schürze



tablier dentaire



delantal dental



歯科エプロン

VISITOR KIT



- Ref no.:1.1.3.29.1

1. Clip Cap: PP10gsm, Size 21" --1 pc
2. Shoe Cover: PP25gsm, 17*38cm, --2pcs
3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
4. Isolation Gown: PP25gsm, 130*140cm --1pc

Packing: 1kit/bag, 20bags/ctn, 29*25*30cm

The components can be customized

			
Besucher-Kit	kit de visiteur	kit de visitante	訪問者キット

VISITOR KIT



- Ref no.:1.1.3.29.2

1. Clip Cap: PP10gsm, Size 21" --1 pc
2. Shoe Cover: PP25gsm, 17*38cm, --2pcs
3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
4. Visitor gown: PP25gsm, 130*140cm --1pc

Packing: 1kit/bag, 20bags/ctn, 29*25*30cm

The components can be customized

			
Besucher-Kit	kit de visiteur	kit de visitante	訪問者キット

VISITOR KIT



- Ref no.:1.1.3.29.3

1. Clip Cap: PP10gsm, Size 21" --1 pc
2. Shoe Cover: PP25gsm, 17*38cm, --2pcs
3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
4. Visitor gown: PE30g/pc, 130*140cm --1pc

Packing: 1kit/bag, 20bags/ctn, 29*25*30cm

The components can be customized

			
Besucher-Kit	kit de visiteur	kit de visitante	訪問者キット

Mother KIT



- Ref no.:1.1.3.29.4

1. Clip Cap: PP10gsm, Size 21" --1 pc
2. Shoe Cover: PP25gsm, 17*38cm, --2pcs
3. Face Mask: PP18/20/25gsm, 17.5*9.5cm --1pc
4. Mother Shirt: Spunlace 50gsm+SMS 35gsm, 130*140cm --1pc

Packing: 1kit/bag, 20bags/ctn, 29*25*30cm

The components can be customized

			
Mutter-Kit	kit mère	kit de la madre	母キット