

Surgical aspiration with cannula



- Ref no:1.4.1
- Description:Suction device composed of a surgical cannula pre-installed
- Specification:Length 2.70m, Internal diameter 0.6cm
- Different size available
- Material:PVC
- Packing: 1pc/blister, 20pcs/box, 40pcs/ctn

			
Absaugkanüle	canule d'aspiration	cánula de aspiración	吸引カニューレ

Surgical aspiration tube



- Ref no:1.4.2
- Description:Transparent PVC suction tube with 2 connectors at the extremities
- Different size available
- Specification:Length 2.70m, Internal diameter 0.6cm
- Material:PVC
- Packing: 1pc/blister, 20pcs/box, 40pcs/ctn

			
Absaugkanüle	canule d'aspiration	cánula de aspiración	吸引カニューレ

Surgical cannula



- Ref no:1.4.3
- Description:Specially for implantologists, it can be fixed to any standard PVC tube
- Specification:Length 25cm, Ø4/6/8mm tip, Crown/Plain Tip,
- With/Without Vent
- Material:PVC
- .-Packing: 1pc/blister, 20pcs/box, 40pcs/ctn

			
Chirurgische Kanüle	canule chirurgicale	cánula quirúrgica	外科カニューレ


Surgical cannula



- Ref. no.:1.4.4
- Description: Precise and effective suction cannula, which reaches the narrowest areas
- Specification: Length 18cm, Ø2.5mm tip
- Material: PVC
- Packing: 1pc/blister, 20pcs/box, 40pcs/ctn

			
Chirurgische Kanüle	canule chirurgicale	cánula quirúrgica	外科カニューレ

Suction Catheter



- Ref no : 1.4.5
- Made of non-toxic medical grade PVC
- Transparent, mist surface is available
- Supplied sterile in individual peelable polybag or blister pack
- Size: F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- Sterile
- Packing: 1pc/bag, 100pcs/ctn or as your required

			
Absaugkatheter	Cathéter d'aspiration	catéter de succión	吸引カテーテル

Stomach Tube



- Ref no : 1.4.6
- Transparent, with X-ray is available
- Supplied sterile in individual peelable polybag or blister pack
- Size: F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- Sterile
- Packing: 1pc/bag, 100pcs/ctn or as your required

			
Magensonde	estomac Tube	tubo de estómago	胃管

Nelaton Catheter



- Ref no : 1.4.7
- Soft, frosted and kink resistant PVC tubing
- Catheters are manufactured from non-toxic, non-irritant, soft PVC compatible with Catheter lubricants
- Colour coded plain connector for easy identification of size
- Sterile product for single use
- Supplied sterile in individual peelable polybag or blister pack
- Size: F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- Packing: 1pc/bag, 100pcs/ctn or as your required

			
Nelaton Katheter	Nélaton Catheter	Nelaton catéter	ネラトンカテーテル

Rectal Tube



- Ref no : 1.4.8
- Inserted into the rectum in order to relieve flatulence which has been chronic and which has not been alleviated by other methods
- Made of Non-toxic PVC, medical grade
- Supplied sterile in individual peelable polybag or blister pack
- Size: F24, F26, F28, F30, F32, F34, F36
- Packing: 1pc/bag, 100pcs/ctn or as your required

			
Darmrohr	Tube rectal	tubo rectal	直腸チューブ

Feeding Tube



- Ref no : 1.4.9
- A safe, easy method for feeding weak patient.
- Made of non-toxic medical grade PVC
- Transparent, with X-ray is available
- Supplied sterile in individual peelable polybag or blister pack
- Size: F4, F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- Packing: 1pc/bag, 100pcs/ctn or as your required

			
Ernährungssonde	Sonde d'alimentation	Tubo de alimentación	栄養チューブ

Y irrigation set



- Ref no : 1.4.10
 - PVC material ,Y irrigation ,
 - tube of 150cm+30cm,2way
 - with flow regulator
- Packing:1pc/bag,100pcs/ctn or as your required

			
Bewässerungsset	irrigation Set	Conjunto de riego	灌漑セット


Mucus Extractor



- Ref no : 1.4.11
 - For aspiration of secretion from oropharynx in newly born babies to ensure free respiration.
 - Made of Non-toxic PVC ,medical grade
 - Supplied sterile in individual peelable polybag or blister pack
 - Size:,F8,F10,F12
 - Sterile
- Packing:1pc/bag,50pcs/ctn or as your required

			
Ernährungs-sonde	Sonde d'alimentation	Tubo de alimentación	栄養チューブ


2-Way Latex Foley Catheter



- Ref no : 1.4.12
- 2-Way Standard
- Size:FR12-FR26
- Made from natural latex, Silicone coated
- Length: 400mm
- Two kinds of valve available with rubber valve and plastic valve
- Packing:1pc/blister bag,100pcs/ctn

			
Foley-Katheter	Cathéter de Foley	Sonda de Foley	フォーリーカテーテル

3-Way Latex Foley Catheter



- Ref no : 1.4.13
- Size: FR162-FR26
- 3-Way Standard
- Made from natural latex, Silicone coated
- Length: 400mm
- Two kinds of valve available with rubber valve and plastic valve
- Packing: 1pc/blister bag, 100pcs/ctn

			
Foley-Katheter	Cathéter de Foley	Sonda de Foley	フォーリーカテーテル

Urethral Catheter (Red Latex)



- Ref no : 1.4.14
- Made from natural latex, Silicone coated
- Length: 400mm
- Two kinds of valve available with rubber valve and plastic valve
- Packing: 1pc/blister bag, 100pcs/ctn

			
UrethraKatheter	Catheter urétral	El catéter uretral	尿道カテーテル

Malecot Catheter



- Ref no : 1.4.15
- Made from natural latex, Silicone coated
- Length: 400mm
- Packing: 1pc/blister bag, 100pcs/ctn

			
Malecot-Katheter	Malecot Catheter	malecot catéter	マレコットカテーテル

2-Way Silicone Foley Catheter...



- Ref no : 1.4.16
- 2-Way Standard
- 100% silicone, medical grade
- Length: 400mm
- Packing: 1pc/blister bag, 100pcs/ctn

			
Foley-Katheter	Cathéter de Foley	Sonda de Foley	フォーリーカテーテル

3-Way Silicone Foley Catheter



- Ref no : 1.4.17
- 3-Way Standard
- 100% silicone, medical grade
- Length: 400mm
- Packing: 1pc/blister bag, 100pcs/ctn

			
Foley-Katheter	Cathéter de Foley	Sonda de Foley	フォーリーカテーテル