

Effective Nebulizer



- Ref no:2.2.1
- Maximum compressor pressure: $\geq 0.15\text{Mpa}$
- Compressor free air flow: $\geq 10\text{L/min}$
- Maximum nebulization rate: $\geq 0.20\text{mL/min}$
- Sound level: $\leq 65\text{dB(A)}$
- Patent design for the clamping groove
- Packing: 1sets/box, 2sets/ctn

			
Vernebler Maschinen	Machine à nébulisateur	nebulizadore	ネブライザー機

Air-compressing nebulizer



- Ref no:2.2.2
- Air-compressing nebulizer
- Maximum compressor pressure: $\geq 0.15\text{Mpa}$
- Compressor free air flow: $\geq 10\text{L/min}$
- Maximum nebulization rate: $\geq 0.20\text{mL/min}$
- Sound level: $\leq 65\text{dB(A)}$
- Packing: 1sets/box, 2sets/ctn

			
Vernebler Maschinen	Machine à nébulisateur	nebulizadore	ネブライザー機

Ultrasonic nebulizer (Single Head)



- Ref no:2.2.3
- Particle diameter $1-5\mu\text{m}$
- Max nebulizing rate $\geq 2\text{ml/min}$
- Medicine cup capacity 50ml
- Setting of Nebulization time: Continuously or from 0 to 60 min for time setting
- Product size :29*17*21CM
- Suggested Applications: Clinical, Nursing Home, Home Care

			
Vernebler Maschinen	Machine à nébulisateur	nebulizadore	ネブライザー機