

Single head Stethoscope



- Ref no : 2.1.1
- The Single head stethoscope is preferred by Cardiac Specialists for monitoring heart sounds.
- Preferred by cardiac specialists.
- For adult ,oxidized aluminum or electroplated zinc
- Chest piece with finger rest for firm application
- 22" Latex free tubing
- Packing:1pc/color box,100pcs/ctn

			
Stethoskop	Et sur Kope S Yo	Y en S Yo s Kope	聴診器

Stainless steel stethoscope



- Ref no : 2.1.2
- Special design of binaural for suiting the ear structure. With uniquely designed binaural thick tubing assembly.
- For adult,deluxe type,with inner spring binaural tubing assembly.
- Special binaural design to suit ear structure
- Diameter of head:47.5mm
- Foam box packing
- Packing:1pc/color box,100pcs/ctn

			
Stethoskop	Et sur Kope S Yo	Y en S Yo s Kope	聴診器

Fetal stethoscope



- Ref no : 2.1.3
- Specially designed chestpiece for pregnancy use
- Both checking the heart beat of fetus and the heart beat of the pregnant woman
- by rotating the chestpiece
- Packing: -Packing:1pc/color box,100pcs/ctn

			
Stethoskop	Et sur Kope S Yo	Y en S Yo s Kope	聴診器

Special dual head stethoscope



- Ref no : 2.1.4
- This is the most popular stethoscope.
- Chrome-plated binaurals and fully rotatable chestpiece
For adult and child.
- Internal spring binaural for professional purpose
- PVC zipper case or box packing.
- Packing:1pc/color box,100pcs/ctn

			
Stethoskop	Et sur Kope S Yo	Y en S Yo s Kope	聴診器

Sethoscope for teaching use



- Ref no : 2.1.5
- Patient-friendly non-chill rim on both the diaphragm and the bell
- Snap-tight soft-sealing eartips assure an excellent acoustical seal
in the ear canal, maximum comfort
- Folds easily without kinking or breaking
- With two binaurals
- Packing:1pc/color box,100pcs/ctn

			
Stethoskop	Et sur Kope S Yo	Y en S Yo s Kope	聴診器

Mercurial sphygmomanometer



- Ref no : 2.1.6
- Measure scope:0-300 mmHg
- Accuracy: 3mm Hg
- sub-division:2mm Hg
- Color available :blue,gray red,green,yellow,black,brown
- Available with nylon cuff or cotton cuff
- With single head stethoscope
- Packing:1pc/box,20pcs/ctn

			
Sphygmomanometer	sphyg bousculait à boire r	p monitor de RESIÓN yr	血圧計

Wrist digital blood pressure monitor



- Ref no : 2.1.7
- Display type:LCD35*45cm
- Memory:7 sets including blood pressure,pulse rate,date,time
- Measurement rage:0-33 mmHG(pressure) 40-200 pulse/min(pulse)
- Accuracy: 3mm Hg(pressure) 5%(Pulse)
- Intelligent Measure:Average function for last three measurement
- Pressure detection:pressure sensor
- Power supply:2 batteries(AAAX2)
- Battery management:3 minutes power-off/low battery indicator
- Product size:72*69*27mm
- Packing:1pc/color box,40pcs/ctn

			
Sphygmomanometer	sphyg bousculait à boire r	p monitor de RESIÓN yr	血圧計

Auto digital blood pressure monitor



- Ref no : 2.1.8
- Display type:LCD52*45cm
- Memory:60 sets including blood pressure,pulse rate,date,time
- Method:oscillometric
- Measurement rage:0-33 mmHG(pressure) 40-200 pulse/min(pulse)
- Accuracy: 3mm Hg(pressure) 5%(Pulse)
- Intelligent Measure:Automatic pressurization,average function for last three measurement
- Depressurization:automatic constant speed exhaust
- Pressure detection:pressure sensor
- Power supply:4 batteries(AAAX4)
- Battery management:3 minutes power-off/low battery indicator
- Product size:135*93*59mm
- Packing:1pc/color box,40pcs/ctn

			
Sphygmomanometer	sphyg bousculait à boire r	p monitor de RESIÓN yr	血圧計


Palm Type Aneroid Sphygmomanometer



- Ref no : 2.1.9
- Measure scope:0-300mmHg
- Accuracy:±3mmHg
- Sub-division:2mmHg
- Nylon cuff with metal ring
- ABS plastic gauge with chrome plated
- Diameter:51mm
- PVC bladder and bulb Leather zipper bag
- 1 pc in a leather bag,then in a gift box,50pcs/ctn

			
Sphygmomanometer	sphyg bousculait à boire r	p monitor de RESIÓN yr	血圧計

Palm Type Aneroid Sphygmomanometer



- Ref no : 2.1.10
- Measure scope:0-300mmHg
- Accuracy:±3mmHg
- Sub-division:2mmHg
- Nylon cuff with silk printing with metal ring
- ABS plastic gauge
- Diameter:69mm PVC bladder and bulb
- Nylon zipper bag
- Color:black,grey,green,red
- 1 pc in a nylon bag,then in a gift box ,50pcs/ctn

			
Sphygmomanometer	sphyg bousculait à boire r	p monitor de RESIÓN yr	血圧計

Pulse Oximeter



- Ref no 2.1.11
- Pulse rate value display, bar graph display
- Low-voltage indication: low-voltage indicator appears before working abnormally which is due to low-voltage
- Automatically power off function: when the device is under the state of measuring interface . it will automatically power off within 5 seconds if the finger falls out of probe
- Various color of cover can be selected

			
Pulsoximeter	Pulse oxymètre	oxímetro de pulso	パルスオキシメータ

Blood Glucose Meter



- Ref no 2.1.12
- Test strip: SXT blood glucose test strip
- Sample: capillary whole blood
- Test range: 2.2mmol/L~27.8mmol/L(40mg/dl~500mg/dl)
- Test time: 25s
- Battery: Two "AAA", 1.5V Alkali battery
- Battery life: about 1000 tests
- Memory: 220, and automatically calculate the average value of day, 14-day and 28-day.
- Test temperature: 10 ~35
- Dimension: 85(L) x 67(W) x 25(H)mm , Weight: about 100g

			
Blutzuckermessgerät	B O o d G U k o s e U m a t e r	Medidor de glucosa	血糖計

Pulse Oximeter



- Ref no 2.1.13.1
- Small in volume、light in weight and convenient in carrying
- Operation of the product is simple, low power consumption
- Operation menu for the function setting
- SpO2 value display
- Pulse rate value display, bar graph display
- Pulse waveform display
- Screen brightness can be changed
- Pulse rate sound indication

			
Pulsoximeter	Pulse oxymètre	oxímetro de pulso	パルスオキシメータ

Penlight



- Ref no : 2.1.13.2
- 100% Brand New. Weight: 20g
- Pen length: 13.5cm ,Bright clear light
- Pupil gauge: 2mm - 9mm
- Convenient push-button on/off switch
- Low power consumption and long use time
- Powered by 2 x 1.5 AAA battery.(Not included)
- Packing: 6pcs/bag, 120pcs/box, 600pcs/ctn

			
Pen l ght	Pen l ght	Pluma que lucha	ペンライト

Pen Light



- Ref no : 2.1.13.3
- light color:yellow,material:ABS,with battery
- Packing:1set/polybag,40set/mid box,480/ctn

			
Pen l ght	Pen l ght	Pluma que lucha	ペンライト

Reflex Hammer



- Ref no : 2.1.14.1
- Multifunctional reflex hammer
- Weight: 65g
- Length: 150mm Width: 45mm
- Stainless steel handle
- Different color is available
- Packing:1pc/box,100pcs/ctn

			
Re f tut x Hamme r	Re Tut f x r Hamme	martillo de reflejos	反射ハンマー

Infrared Ear Thermometer



- Ref no : 2.1.15
- Display range: 34.0~44.0oC/93.2~111.2oC
- Accuracy: $\pm 0.2^{\circ}\text{C}/\pm 0.4^{\circ}\text{F}$ Min. scale: 0.1
- Memory: 10 memories
- Measurement time (reference only): 1 second
- Auto shut-off (1min after idle)
- Battery: AAA x 2, alkaline batteries
- Battery life: approx. 4,000 takes
- Size: 134 x 70 x 34mm
- Requires new disposable probe cover after each use
- Packing:1pc/box,20pcs/ctn

			
Thermometer r	thermomètre	termómetro	温度計

Digital Thermometer



- Ref no : 2.1.16
- Display range: 32.0~43.0°C
- Accuracy: $\pm 0.1^{\circ}\text{C}$
- Min. scale: 0.1
- Measurement time (reference only):
 - a) 60 ± 10 seconds (oral/rectal)
 - b) 100 ± 20 seconds (underarm)
- Beeper function Auto shut-off
- Battery: 1.5V button battery (LR/SR-41)
- LCD:15.5*6.5mm
- 1pc/box,200pcs/ctn

			
Thermometer r	thermomètre	termómetro	温度計

Oral Thermometer



- Ref no : 2.1.17
- For oral use
- With plastic case packing
- Single scale 32.1-42.0
- Double scale 32.0-42.0(89.6-107.6)
- Accuracy: 0.1(0.2)
- Taking time:3 minutes
- Packing:12pcs/box,720pcs/ctn

			
Thermometer r	thermomètre	termómetro	温度計

Armpit Thermometer



- Ref no : 2.1.18
- For armpit use
- With deluxe plastic case packing
- Packing:12pcs/box,720pcs/ctn

			
Thermometer r	thermomètre	termómetro	温度計