

Disposable Electrode



- Ref. no. 1.16.1
- Liquid glue or Hardened Hydrogel available
- Non-woven Fabrics backing or PE foam available
- Size: $\Phi 55\text{mm}$
- Form: circular
- Adult
- Packing: 50pcs/bag, 2000pcs/ctn

			
Einweg-Elektrode	Electrode jetable	En Dji s Posa b tut Sí filtrado ct	ディスボ電極

Disposable Electrode



- Ref. no. 1.16.2
- Suitable for normal ECG tests in diagnostic or monitoring, it uses an Ag/AgCl sensor element and solid conductive & adhesive Hydro-gel or wet conductive gel for adhesion.
- Liquid glue or Hardened Hydrogel available
- Non-woven Fabrics backing or PE foam available
- Size: 43*49mm(Adult/Child), 36*42mm(Child), 30*36mm(infant)
- Form: drop shaped
- Packing: 50pcs/bag, 2000pcs/ctn

			
Einweg-Elektrode	Electrode jetable	En Dji s Posa b tut Sí filtrado ct	ディスボ電極

Disposable Electrode



- Ref. no. 1.16.3
- Suitable for normal ECG tests in diagnostic or monitoring, it uses an Ag/AgCl sensor element and solid conductive & adhesive
- Hydro-gel or wet conductive gel for adhesion.
- Silver-Silver chloride, self adhesive, PE foam backing
- Snap fastener, 2.8 x 4.5 cm,
- Packing: 50pcs/bag, 2000pcs/ctn

			
Einweg-Elektrode	Electrode jetable	En Dji s Posa b tut Sí filtrado ct	ディスボ電極


Disposable Electrosurgical Pencil



- Ref. no. 1.16.4
- Disposable Electrosurgical Pencil(ESU Pencil) is for general surgery, chest surgery usually.
- High-quality ABS shaped
- 100%% EO gas sterilized packing
- plug standard connector
- 3 meters cable, 2.5m and 5m is for slectable
- Electrode Exposed Length :40mm
- Function: Cut and Coagulation
- Packing:1pc/blister,100pcs/ctn

			
Einweg-Elektrode	Electrode jetable	En Dji s Posa b tut Sí filtrado ct	デイスボ電極

ECG paper



- Ref no : 1.16.5
- Different medical paper for EKG /ECG/EEG/CTG/USG checking in Hospital
- Various size available.

			
EKG-Papier	papier ECG	papel del ECG	ECG 紙