

Non woven face Mask-2ply



- Ref no. : 1.1.5.1.1.1
- Material: Non woven
- 2-ply (pp non woven+ pp non woven) , with tie on
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Non woven face Mask -2ply



- Ref no. : 1.1.5.1.1.2
- Material: Non woven
- 2-ply (pp non woven+ pp non woven) , with elastic earloop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Non woven face Mask-3ply



- Ref no. : 1.1.5.1.2.1
- Material: Non woven
- 3-ply (pp non woven+ melt blown +pp non woven) , with tie on
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Non woven face Mask-3ply



-Ref no. : 1.1.5.1.2.2

Material: Non woven

3-ply (pp non woven+ melt blown +pp non woven) , with elastic ear loop

-Size: 17.5*9.5cm

-Packing: 50pcs/box, 40boxes/ctn, packing can be customized

- Color option: white, color could be customized



non-woven facemask



non tissé masque



no tejida máscara facial



不織布フェイスマスク

Non woven face Mask-4ply



-Ref no. : 1.1.5.1.3.1

- Material: Non woven

- 4-ply (pp non woven+ pp non woven+ melt blown +pp non woven) , with tie on

-Size: 17.5*9.5cm

-Packing: 50pcs/box, 40boxes/ctn, packing can be customized

- Color option: white, color could be customized



non-woven facemask



non tissé masque



no tejida máscara facial



不織布フェイスマスク

Active Carbon Face Mask



-Ref no. : 1.1.5.1.3.2

- Material: SBPP

- With tie on

-Size: 17.5*9.5cm

-Packing: 50pcs/box, 40boxes/ctn, packing can be customized

- Color option: Black, color could be customized



Aktivkohle facemask



un masque facial de charbon
actif



mascarilla de carbón activo



活性炭フェイスマスク

Active Carbon Face Mask



- Ref no. : 1.1.5.1.3.3
- Material: SBPP
- With elastic ear loop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: black, color could be customized

			
Aktivkohle facemask	un masque facial de charbon actif	mascarilla de carbón activo	活性炭フェイスマスク

Non woven face Mask-4ply



- Ref no. : 1.1.5.2.1
- Material: Non woven
- 4-ply (pp non woven+ pp non woven+ melt blown +pp non woven) , with elastic earloop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Type II face Mask



- Ref no. : 1.1.5.3
- Material: Non woven
- 3-ply (pp non woven+ melt blown +pp non woven) , with elastic ear loop, 99% bacterial filtration efficiency
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Type II R face Mask



- Ref no. : 1.1.5.4
- Material: Non woven
- 3-ply (pp non woven+ melt blown +pp non woven) , with elastic ear loop, 99% bacterial filtration efficiency + 120 of synthetic blood penetration resistance
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク


Cellulose face Mask



- Ref no. : 1.1.5.5
- Material: Cellulose +Non woven
- 3-ply (Cellulose+ melt blown +pp non woven) , with elastic ear loop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: Yellow, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Wetland face Mask



- Ref no. : 1.1.5.6
- Material: Wetland paper +Non woven
- 3-ply (Wetland paper+ melt blown +pp non woven) , with elastic ear loop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 40boxes/ctn, packing can be customized
- Color option: yellow, color could be customized

			
non-woven facemask	non tissé masque	no tejida máscara facial	不織布フェイスマスク

Face Mask with shield



- Ref no.:1.1.5.7.1
- Material: Non woven
- With shield, with tie on
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 10boxes/ctn, packing can be customized
- Color option: white or blue color can be customized

			
facemask mit Schild	masque avec le bouclier	mascarilla con escudo	シールド付きフェイスマスク

Face Mask with shield



- Ref no.:1.1.5.7.2
- Material: Non woven
- With shield, with elastic ear loop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 10boxes/ctn, packing can be customized
- Color option: white or blue color can be customized

			
facemask mit Schild	masque avec le bouclier	mascarilla con escudo	シールド付きフェイスマスク

Paper face Mask-1ply



- Ref no. : 1.1.5.8.1
- Material: Paper
- 1-ply paper
- Size: 20*7cm
- Packing: 100pcs/box, 100boxes/ctn, packing can be customized
- Color option: white

			
Papier facemask	papier masque	mascarilla de papel	紙のフェイスマスク

Paper face Mask-2ply



- Ref no. : 1.1.5.8.2
- Material: Paper
- 2-ply paper
- Size: 20*7cm
- Packing: 100pcs/box, 100boxes/ctn, packing can be customized
- Color option: white

			
Papier facemask	papier masque	maskarilla de papel	紙のフェイスマスク

Protective FFP1 face Mask



- Ref no. : 1.1.5.9.1
- Material: Non woven
- FFP1 standard, cone type and without valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP1 face Mask



- Ref no. : 1.1.5.9.2
- Material: Non woven
- FFP1 standard, cone type and with valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP1 face Mask



- Ref no. : 1.1.5.9.3
- Material: Non woven
- FFP1 standard, boat type and without valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP1 face Mask



- Ref no. : 1.1.5.9.4
- Material: Non woven
- FFP1 standard, boat type and with valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP2 face Mask



- Ref no. : 1.1.5.10.1
- Material: Non woven
- FFP2 standard, cone type and without valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP2 face Mask



- Ref no. : 1.1.5.10.2
- Material: Non woven
- FFP2 standard, cone type and with valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP2 face Mask



- Ref no. : 1.1.5.10.3
- Material: Non woven
- FFP2 standard, boat type and without valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP2 face Mask



- Ref no. : 1.1.5.10.4
- Material: Non woven
- FFP2 standard, boat type and with valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP3 face Mask



- Ref no. : 1.1.5.11.1
- Material: Non woven
- FFP3 standard, cone type and without valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP3 face Mask



- Ref no. : 1.1.5.11.2
- Material: Non woven
- FFP3 standard, cone type and with valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP3 face Mask



- Ref no. : 1.1.5.11.3
- Material: Non woven
- FFP3 standard, boat type and without valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective FFP3 face Mask



- Ref no. : 1.1.5.11.4
- Material: Non woven
FFP3 standard, boat type and with valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective N95 face Mask



- Ref no. : 1.1.5.12.1
- Material: Non woven
N95 standard, cone type and without valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective N95 face Mask



- Ref no. : 1.1.5.12.2
- Material: Non woven
N95 standard, cone type and with valve
- Size: 13*13cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective N95 face Mask



- Ref no. : 1.1.5.12.3
- Material: Non woven
- N95 standard, boat type and without valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Protective N95 face Mask



- Ref no. : 1.1.5.12.4
- Material: Non woven
- N95 standard, boat type and with valve
- Size: 20*10cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Molded Mask



- Ref no. : 1.1.5.13
- Material: Non woven
- Cone type
- Size: 11.5*12cm
- Packing: 20pcs/box, 20boxes/ctn, packing can be customized
- Color option: white, color could be customized

			
Schutz facemask	masque de protection	no tejida máscara facial	保護フェイスマスク

Clear mask



- Ref no. : 1.1.5.14
- Material: PC
- Size: 140*60mm
- Packing: 10pcs/box, 50boxes/ctn, packing can be customized
- Color option: transparent

			
C Ea r oder sk	C Ea r ou sk	C Ea r o sk	クリアマスク

Anti-fog Face Shield



- Ref no.:1.1.5.15
- Material: PET with elastic band
- Size:22*33cm
- Packing: 10pcs/Plastic bag, 20bags/ctn, can be customized

			
Anti-Fog-Face Shield	Anti-brouillard Face Shield	Antivaho careta	アンチフォグフェイスシールド

Face Mask with shield,with foam



- Ref no.:1.1.5.16
- Material: Non woven
- With shield, foam, earloop
- Size: 17.5*9.5cm
- Packing: 50pcs/box, 10boxes/ctn, packing can be customized
- Color option: white/yellow/blue, color can be customized

			
facemask mit Schild	masque avec le bouclier	mascarilla con escudo	シールド付きフェイスマスク

Beard cover



- Ref no. : 1.1.5.17
- Material: Non woven
- with single elastic
- Size: 20*40cm, could be customized
- Packing: 100pcs/bag, 10bags/ctn, packing can be customized
- Color option: white, color could be customized



Bart Abdeckung



couverture barbe



cubierta de la barba



ひげカバー

Beard cover



- Ref no. : 1.1.5.18
- Material: Non woven
- with double elastic
- Size: 20*40cm, could be customized
- Packing: 100pcs/bag, 10bags/ctn, packing can be customized
- Color option: white, color could be customized



Bart Abdeckung



couverture barbe



cubierta de la barba



ひげカバー