

Stainless Steel Bedside Cabinet



- Ref.no.4.2.1
- Item is made of stainless steel plate;
- One sliding table, one drawer and one locker with galvanized handle;
- Towel hangers on both sides;
- With a removable thermos bottle rack,
which is made of stainless steel rectangular tube;
- Size:L480*W410*H800mm

			
Base-d Herr Kabine t	Base de d M. Kabine t	Base d Sr. Kabine t	ベッドサイドキャビネット

Steel Bedside Cabinet(ABS Top Panel)



- Ref.no. 4.2.2
- Steel Bedside Cabinet(ABS Top Panel)
- Material: Epoxy powder coated steel body and ABS surface
- With one dining board, one drawer and two cabinets(one big and one small)
- With 2 folding racks for hanging towel or other things
- Size: L450*440*900mm

			
Base-d Herr Kabine t	Base de d M. Kabine t	Base d Sr. Kabine t	ベッドサイドキャビネット

Steel Bedside Cabinet(ABS Top Panel)



- Ref.no. 4.2.3
- Steel Bedside Cabinet (ABS Top Panel)
- Material: Epoxy powder coated steel body and ABS surface
- With one dining board, one drawer and two cabinets(one big and one small)
- With 2 folding racks for hanging towel or other things
- Size: L450*440*900mm

			
Base-d Herr Kabine t	Base de d M. Kabine t	Base d Sr. Kabine t	ベッドサイドキャビネット

Luxury Medicine Trolley



- Ref.no.4.2.4
- Size:L770*W550*H1000mm

			
Medizin Trolley	Medicine Trolley	Mejishine T y Ro'e	医学トロリー