

High pressure injector-single head



- Ref. no.4.6.1
- Flow rate: 0.1-10mL/s, increments 0.1mL/s
- Pressure limits: 100-250psi, increments 50psi
- KVO: 0.1mL/s
- Scan delay: 0-600s
- Injection delay: 0-600s
- Max. 8 phases support different injection combination
- 200mL large volume syringe meets most application needs



Hallo gh D ruck Ince c und r



Salut gh p resion Ince c et r



Hola gh p resión Ince c y r



高圧インジェクタ

Dental Surgery Microscope



- Ref. no 4.6.2
- Main double ocular
- The Big objective lens(with micro-adjusting)
- Five steps drum wheels are changed into 3.4x, 5.4x, 8.5x, 13.5x, 21.2x , Floor stand
- The good incline system with damping The good up-down system with spring arm



Den wurde l suchen c Filtration
s Kope



Microscopio dental



Den estaba buscando l c
filtración s Kope



歯科顕微鏡

Vision Video Laryngoscope



- Ref. no .4.6.3
- Fully Portable Video laryngoscope, 210g
- Anti-fog
- Can be connected to bigger monitor and computer
- Real-time take photo and video
- Reusable 304 stainless steel blade and disposable plastic MAC 3 blade!



Koralle s Kope



Kope de corail



coral ry s Kope



喉頭鏡

Dental Lab Induction Heated Casting Machine



- Ref. no .4.6.4
- Middle Frequency Induction Casting Machine is upgraded from Power Supply: 220V/110V 50/60Hz
- Power Consumption: 2,500W, Motor: 370W
- Centrifugal Rev: 550rpm
- Centrifugal Radius: 21cm
- Casting Method: Centrifugal Casting
- Max. Melting Amount: 50g
- Heating Method: Medium Frequency Induction Heating



Den war l oder s Laterne g Matinee



Den était l ou s lanterne g
matinee



máquina de colada dental



歯科鑄造機

Examination Light



- Ref. no .4.6.5
- Voltage:AC90-240V
- Power:LED21W
- Illumination:80000lx
- Weight:5kgs
- Life:50000h



Untersuchungsleuchte



Examen Lumière



examen Luz



検査ライト

36 Holes LED Ceiling Examination Light



- Ref. no .4.6.6
- Voltage:AC90-240V
- Power:LED108W
- Installation height:2800mm-3100mm
- Life:50000h
- Weight:35kgs



Untersuchungsleuchte



Examen Lumière



examen Luz



検査ライト

High Frequency Digital X ray Diagnostic System



- Ref. no .4.6.7
- 50kw high frequency technology,output constant high voltage direct current,low ripple factor,clearer and sharper image
- Automatic Program radiography,provide thousand of type radiographic conditions,can set and save radiographic conditions for the especial requirements.\

			
X ray-Diagnosesystem	X Système de diagnostic ray	Sistema de diagnóstico rayos X	X 線診断システム

Medical Bucky stand



- Ref. no .4.6.8
- Ascending and Descending Range: 560-1625mm
- Film Size: 5"×7" - 17"×17"
- With Imported Korean JPI Grid
- Grid Size: 5"×7" - 18"×18"
- Grid Ratio: 10:1
- Grid Density: 103L/inch
- Focal Distance: 150cm

			
Medical Bucky stehen	Médical Bucky debout	Médico Bucky reposar	医療バックースタンド