

Blood collection tube



- Ref no : 1.12.1
- Specification: 1-10ml
- Material: PET or GLASS
- Packing: 250pcs/bag, 10bag/ctn
- Color option: Transparent
- Different tube available

			
Blutentnahmeröhrchen	Tube de prélèvement de sang	tubocolelector de sangre	採血管

Test tubes



- Ref no : 1.12.3
- Material :PS
- Specification:14*41mm,14*49mm,14*55mm,14*53mm
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Reagenzgläser	Des tubes à essai	Tubos de ensayo	試験管


Pipette



- Ref no : 1.12.2.7
- Specification:3ml
- Material :LDPE
- Packing:10000pcs/ctn
- Color option: Transparent

			
Pipette	Pipette	Pipeta	ピペット

Microscope slide



- Ref no : 1.12.6
- Specification: 76.2 MM. X 25.4 MM. (3" X 1").
- Thickness: 1 mm. - 1.2 mm.
- Packing: 72pcs/box, 2,880pcs/ctn
- Color option: transparent.

			
Objekträger	Microscope slide	portaobjetos de microscopio	顕微鏡スライド

Urine tube with screw cap



- Ref no : 1.12.5
- Specification: 12ml, With a watertight closure
- Material: PS
- Packing: 200pcs/bag, 10bags/ctn
- Color option: transparent

			
Urinröhrchen mit Schraubverschluss	tube d'urine avec bouchon à vis	tubo de la orina con tapón de rosca	スクリーキャップと尿管

Cryovials tubes



- Ref. no.: 1.12.4.14
- 4.5 ml
- Material: polycarbonate
- Specification: conical internal base, o-ring cap.
- Packing: 500pcs/bag, 8bags/ctn
- Color of cap: various for option

			
Rohre Kryovials	tubes cryovials	tubos crioviales	チューブクライオバイアル

Cryovials tubes



- Ref. no.: 1.12.4.13
- 1.5 ml
- Material: polycarbonate
- Specification: conical internal base, o-ring cap.
- Packing: 1000 pcs/bag, 10 bags/ctn
- Color of cap: various for option

			
Rohre Kryovials	tubes cryovials	tubos crioviales	チューブクライオバイアル

Centrifuge Tube



- Ref. no.: 1.12.1.10
- 10ml Centrifuge Tubes
- Conical base
- Packing: 250 pcs/box, 20 boxes/ctn
- Color option: transparent

			
Dazu t Rifuge Tsube	Pour ce faire t Rifuge Tsube	tubo de centrifuga	遠心分離管

PCR tube



- Ref. no.: 1.12.5
- Size: 0.2ml thin wall
- Material: PS
- Packing: 1000 pcs/bag, 10 bags/ctn
- Color option: transparent

			
PCR Tsube	Tsube PCR	tubo de PCR	PCR チューブ

Deep pore plate



- Ref. no.: 1.12.16
- Material: PS
- Size : 1.3ml 96-hole deep pore plate ,
- Specification: Round hole ,U-bottom,Sterilization
- Packing:5pcs/bag,10bags/ctn
- Color option: transparent

			
Porentiefe Platte	Plaque de pores en profondeur	Placa profunda de los poros	ディープポアプレート

Petri Dish



- Ref no : 1.12.2.1
- Material :High grade PS
- Specification: 60ml,90ml
- Packing:1pcs/Plastic bag,400bags/ctn
- Color option: Transparent

			
Petrishale	Boîte de Pétri	Placa de Petri	ペトリ皿

Cuvette cup



- Ref no : 1.12.4.1
- Material: High grade PS
- Multi-sizes for choice
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キューベットカップ

Cuvette cup



- Ref no : 1.12.4.3
- Material :High grade PS
- Size: Multi-sizes for choice
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.4
- Material :High grade PS
- Specification:TECO Blood Clot Meter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.5
- Material :High grade PS
- Systems CA530,CA6000,CA100
- Specification: TECO Blood Clot Meter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.6
- Material :High grade PS
- Specification Systems for Spanish, German Labor
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.7
- Material: High grade PS
- Specification: Trolley for Pacific and German TECO 4
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.8
- Material :High grade PS
- French Organon, zhongshanpeikang
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.9
- Material :High grade PS
- For American Coulter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.10
- Material :High grade PS
- Specification ForTianjinPulisheng and French STAGO
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.11
- Material :High grade PS
- Specification For beckman
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.12
- Material :High grade PS
- Specification: for Beijing Periong Swiss Serozyme
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.13
- Material :High grade PS
- Specification for Roche
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ


Cuvette cup



- Ref no : 1.12.4.14
- Material :High grade PS
- Specification: for Roche, special protein analyzer
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.15
- Material :High grade PS
- Specification for French Biomerieux
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Cuvette cup



- Ref no : 1.12.4.16
- Material :High grade PS
- Specification for American Abbott
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cuvette Tasse	cuvetasse	tazacubeta	キュベットカップ

Biochemical Cup



- Ref no : 1.12.3.1
- Material :High grade PS
- Specification for Japan Dongya F/CIS series Ø30×70
- Japan Kohden/Botestseries Ø30 60
- Mindray BC-2000 Ø30×52
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Biochemische Cup	Coupe Biochemical	Copa bioquímica	生化学カップ

Biochemical Cup



- Ref no : 1.12.3.2
- Material :High grade PS
- Specification for Hitachi 7020,7150 biochemical apparatus
- Beckman CS series
- 700 US Beckman ,Dubang
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Biochemische Cup	Coupe Biochemical	Copa bioquímica	生化学カップ

Biochemical Cup



- Ref no : 1.12.3.3
- Material :High grade PS
- Specification: for Hitachi 7020,7150 biochemical apparatus
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Biochemische Cup	Coupe Biochemical	Copa bioquímica	生化学カップ

Biochemical Cup



- Ref no : 1.12.3.4
- Material :High grade PP or ABS
- Specification :for swiss rose biochemical apparatus
- For Aulinbas AU Reaction cup
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Biochemische Cup	Coupe Biochemical	Copa bioquímica	生化学カップ

Biochemical Cup



- Ref no : 1.12.1.10.1
- Material :High grade PP
- Specification: action cup
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Biochemische Cup	Coupe Biochemical	Copa bioquímica	生化学カップ

Centrifuge Tube



- Ref no : 1.12.1.10.3
- Material :High grade PP
- Specification Microcentrifuge tube
- Size:1.5ml
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent (orange,red,yellow,blue)

			
Zentrifugenröhrchen	Centrifugeuse Tube	Tubocentrifugo	遠心分離管

Centrifuge Tube



- Ref no : 1.12.1.10.4
- Material :High grade PP
- Specification: Microcentrifuge tube
- Size: 0.5ml yellow,1.5ml red ,2ml blue
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Zentrifugenröhrchen	Centrifugeuse Tube	Tubocentrifugo	遠心分離管

Cryovial Tube



- Ref no : 1.12.1.12
- Material :High grade PP
- Specification: Cryovial Tube
- Size:5ml blue,1.8ml yellow,1ml red
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Cryovial Rohr	Tube cryovial	tubocriovial	クライオバイアルチューブ

Racks for Cryovial Tubes



- Ref no : 1.12.1.14.1
- Material :High grade PP
- Specification: 36 Weel (5ml)
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Racks fürCryovial Tubes	Racks fürCryovial Tubes	Bastidores para criovialTubos	クライオバイアルチューブ用ラック

Racks for Cryovial Tubes



- Ref no : 1.12.1.14.2
- Material :High grade PP
- 25 position
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Racks fürCryovial Tubes	Racks fürCryovial Tubes	Bastidores para criovialTubos	クライオバイアルチューブ用ラック

Racks for Cryovial Tube



- Ref no : 1.12.1.14.3
- Material :High grade PP
- Specification: 50 position
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Racks fürCryovial Tubes	Racks fürCryovial Tubes	Bastidores para criovialTubos	クライオバイアルチューブ用ラック

Racks for Cryovial Tube



- Ref no : 1.12.1.14.4
- Material :High grade PP or PS
- Specification :100 position
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Racks fürCryovial Tubes	Racks fürCryovial Tubes	Bastidores para criovialTubos	クライオバイアルチューブ用ラック


PCR tube



- Ref no : 1.12.1.9.1
- Material :High grade PP
- Specification: PCR Tube 8 strips
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
PCR-Röhrchen	tube PCR	tubo de PCR	PCR チューブ

PCR tube



- Ref no : 1.12.1.9.2
- Material :High grade PS
- Specification: PCR Tube 12 strips
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
PCR-Röhrchen	tube PCR	tubo de PCR	PCR チューブ

PCR Plate



- Ref no : 1.12.1.9.2
- Material :High grade PP
- Specification: PCR plate 96 wells
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
PCR Plate	plaque PCR	Placa PCR	PCR プレート


Pipette tips



- Ref no : 1.12.13
- Material :High grade PS
- Specification: :5-10ul,200ul,300ul,1000ul,5000ul,
- Material: Medical grade PP
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Pipettenspitzen	Pointes de pipettes	puntas de pipeta	ピペットチップ

Tip coupled for rack



- Ref no : 1.12.12.2
- Material :High grade PS
- Specification: Boxed 200ul tip with filter
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Tipp gekoppelt für Rack	Astuce couplé pour crémaillère	Tip acoplada para estante	ヒントラックのために結合されま した


Slot for reagent



- Ref no : 1.12.12.5
- Material :High grade PS
- Specification: 25ml
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Slot für Reagenzien	Emplacement pour réactif	Ranura para el reactivo	試薬用スロット

Rack for tips



- Ref no : 1.12.12.2
- Material :High grade PS
- Specification: 1000 up tip box
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Gestell für Tipps	Rack pour obtenir des conseils	Rack para consejos	ヒントについては、ラック


Test tubes



- Ref no : 1.12.3
- Material :High grade PP
- PS material
- Specification:14*41mm,,14*49mm,14*55mm,14*53mm
- Packing:250pcs/bag ,10bag/ctn
- Color option: Transparent

			
Reagenzgläser	Des tubes à essai	Tubos de ensayo	試験管

Urine cup



- Ref no:1.12.1
- Specification:60ml,90ml,120ml
- Materials: cup is made of PP, lid is made of HDPE
- Non-sterile or sterile
- Clear graduation on the cups
- Transparent cup and different color of the lid such as red,
- with label or without label on the cup
- Screw lid
- Packing:250pcs/ctn

			
Urinbecher	urine tasse	taza de la orina	尿カップ

Stool container



- Ref no:1.12.2
- Specification:30ml,60ml
- Materials: cup is made of PP, lid is made of HDPE
- Clear graduation on the cups
- With spoon
- Non-sterile or sterile
- Screw lid
- Packing:1000pcs/ctn

			
Poo Tasse	Poo tasse	taza de Poo	プーカップ

Medicine Cup Without Lid



- Ref no: 1.12.3.1
- Specifications: 1oz. (30ml) ,
- Materials: PP
- Non-sterile
- Clear graduation on the cup
- Transparent

			
Kleine Medizincup	Petite tasse de médecine	Pequeña taza de la medicina	小さな医学カップ

Medicine Cup Without Lid



- Ref no:1.12.3.2
- Specifications: 2oz.60ml,
- Materials: PP
- Non-sterile
- Clear graduation on the cups
- Transparent

			
Kleine Medizincup	Petite tasse de médecine	Pequeña taza de la medicina	小さな医学カップ

Triangular Measuring Cup



- Ref no:1.12.4.1
- PS, material,
- Capacity: 1000ml
- With graduated marks
- Color:transparent

			
Triangle Tasse	Triangle tasse	la taza del triángulo	三角形の計量カップ

Triangular Measuring Cup



- Ref no:1.12.4.2
- PP, material,
- Capacity: 1000ml
- With graduated marks
- Color:white

			
Triangle Tasse	Triangle tasse	la taza del triángulo	三角形の計量カップ

Test tube brush



- Ref no:1.12.5
- Material: Nylon
- Specification: length 20cm, diameter:2,3,4,5,6,7mm
- Packing: 1pc/bag, 2000pcs/ctn
- Color: white or black

			
Reagenzglasbürste	Brosse pour tubes à essais	Cepillo para tubos de ensayo	試験管ブラシ