

Oxygen regulator



- Ref no:2.4.1.1
- Flow regulation: 1`10L/n
- Buck scope: 2.0`3.0kg
- Safety valve pressure rang: 3.5kg±0.5kg
- Packing:1sets/box,10sets/ctn

			
Kontakt xy Gen Bein und r	Contactez-xy Gen jambe	Contacto xy Gen	酸素調節器

Oxygen regulator



- Ref no:2.4.1.2
- Flow adjustable 1-10L/n ,Input pressure : 0.1-0.12Mpa
- Pressure safety valve : 0.17Mpa
- Wave of devices attached
- Packing:1set/box/2sets/ctn

			
Kontakt xy Gen Bein und r	Contactez-xy Gen jambe	Contacto xy Gen	酸素調節器

Oxygen regulator



- Ref. no.:2.4.1.3
- Input:2200psi,output:50psi-22psi,
- flux:12 levels in adjustment
- Packing:1pc/box,20pcs/ctn

			
Kontakt xy Gen Bein und r	Contactez-xy Gen jambe	Contacto xy Gen	酸素調節器

Portable Phlegm Suction Unit



- Ref no:2.4.2
- High negative pressure, large flux, low noise
- Voltage: AC220V±10%.50Hz±2%
- Power:90VA
- Max negative pressure & #118alue: ≥0.075MPa
- PUmping rate: ≥18L/min
- Reservoir volume:1000mL
- Packing:1sets/box,2sets/ctn

			
Saug-Maschine	Machine d'aspiration	máquina de succión	吸引機

Elisa microplate reader



- Ref no.: 2.4.3
- .Touch screen, Large LCD display
- .Multiple plate formats capability
- .Extensive on-board data analysis
- .Choice of reading speed
- .Microplate shaking
- 6.10000 test results

			
Mikroplatten -Reader	Il a l'air c filtrat p § est rare d	Parece c filtrado p § Es raro	マイクロプレートリーダー
		d	

Portable Ultrasound Digital Scanner



- Ref no.: 2.4.4
- Hospitals and clinics for diagnosing abdomen, urology, Cardiology, OB/GYN, small parts examination and orthopaedics.
- Scanning mode : Convex/Linear/Micro-convex
- Display modes: B,B + B,B + M, 4B
- Scanning Depth: 240mm
- Gray scales: 256
- Cine-loop: 256 frames

			
digitale Scanner	Jigita I S Kanye r	Jigita I S r Kanye	デジタルスキャナ

Foetal Monitor Cardiotacograph



- Ref no.: 2.4.5
- Accurate FHR detection
- High sensitivity (1MHz)
- Automatic foetal movement detection
- A4 size wide chart paper
- Uses inexpensive fax paper
- Intuitive user interface
- Easy to upgrade

			
Fetale-Monitor	moniteur foetal	Monitor fetal	胎児モニタ

Dental Ultrasonic Cleaner



- Ref no.:2.4.6.1
- Dentures: Retainers, False Teeth, Invisalign Braces, outhguards, etc. Jewelry: Earrings, Rings, Bracelets, Coins, etc.
- Ultrasonic Frequency 42,000 Hz
- Tank Material Stainless Steel
- Tank Capacity 1 Cup (140 mL): 3 x 3 x 1 inch Tank Size
- Timer 5 minute cleaning cycle
- Power Supply AC 100 ~ 120V, 50 / 60Hz
- Power 20 W
- Unit Size 7.12 X 4.8 X 4.8 inches (L x W x H) ,N.W. 1lb 7oz

			
Zahnreiniger	nettoyant dentaire	limpiador dental	歯科クリーナー

Dental Ultrasonic Cleaner



- Ref no.:2.4.6.2
- Ultrasonic Frequency 42,000 Hz
- Tank Material Stainless Steel
- Tank Capacit:600ml
- Timer 5 minute cleaning cycle
- Power Supply :AC220-240V,50/60HZ
- or AC100-120V,50/60HZ
- Inner size : 155*95*52MM
- Ourter size : 206*152*132MM Weight:1.35kgs

			
Zahnreiniger	nettoyant dentaire	limpiador dental	歯科クリーナー

Dental Autoclave



- Ref no.:2.4.7
- No vacuum, sterilizing the solid instrument.
- LCD display,2 programs for your choose 3.Full-Automatic operation and door lock system, as well as the feeding water system.
- can show the reason of failure on LCD, do not need to check recording to code menu

			
Dental Autoklav	dentaire Autoclave	dental Autoclave	歯科オートクレーブ

Difibrillator



- Ref no.:2.4.8
- Energy : 0, 3, 5, 7, 10, 20,30, 50,100,200,300,360 joule
- Charging time : At 360J less than 10 sec
- Paddle options : Reusable external adult and pediatric paddles
- Operates from the mains (AC) and with its internal Lead-acid battery

			
Defibrillator	défibrillateur	desfibrilador	除細動器