


Nasal cannula-Soft prong



- Ref no:1.14.1.1
- provide flexibility and a more customized fit for patients
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm and 6mm
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

			
Nasenkanüle	canule nasale	cánula nasal	鼻カニューレ

Nasal cannula-Straight prong



- Ref no:1.14.1.2
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm or 6mm
- Soft tubing to increase patient comfort
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

			
Nasenkanüle	canule nasale	cánula nasal	鼻カニューレ

Nasal cannula-For mouth and nose



- Ref no:1.14.1.3
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm or 6mm
- provide flexibility and a more customized fit for patients
- Soft tubing to increase patient comfort
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

			
Nasenkanüle	canule nasale	cánula nasal	鼻カニューレ

Nasal cannula-For mouth and nose



- Ref no:1.14.1.4
- Color:Clear or green
- Matrial:soft PVC and regular PVC
- Diameter:5mm or 6mm
- provide flexibility and a more customized fit for patients
- Soft tubing to increase patient comfort
- Provides crush-proof and kink-proof tubing for continuous flow
- Packing:1pc/polybag,50pcs/ctn or as your required

			
Nasenkanüle	canule nasale	cánula nasal	鼻カニューレ

Nasal oxygen catheter



- Ref no: 1.14.1.5
- Size:FR8,fr10,fr12,fr14
- Suitable for easy and efficient administration of oxygen.
- Sof and smooth open distal end facilitates non-traumatic insertion.
- Funnel connector at proximal end.
- Three type are available,with moveable sponge,with non-moveable
- sponge at end,and without sponge.

			
Nasenkanüle	canule nasale	cánula nasal	鼻カニューレ

O2+CO2 Sampling Cannula



- Ref no: 1.14.1.6
- Adult or Child prong
- Color:Clear or green
- Diameter:5mm or 6mm
- Designed for Oxygen supply and Sampling of exhaled co2 gas simultaneously
- Soft prongs minimize irritation

			
Sampling Cannula	échantillonnage Cannula	La cánula de muestreo	サンプリングカニューレ

Nebulizer mask



- Ref no: 1.14.2.1
- 6ML capacity anti-spill jar to prevent leakage
- Color: Clear or green
- Diameter: 5mm or 6mm
- Kink-resistant
- Elastic strap fixation
- Latex-free
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Vernebler Maske	masque nébuliseur	máscara nebulizador	ネブライザーマスク

Nebulizer mask



- Ref no: 1.14.2.2
- 15ML capacity anti-spill jar to prevent leakage
- Color: Clear or green
- Diameter: 5mm or 6mm
- Kink-resistant
- Elastic strap fixation
- Latex-free
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Vernebler Maske	masque nébuliseur	máscara nebulizador	ネブライザーマスク

Nebulizer with mouth piece



- Ref no: 1.14.2.3
- 15ML or 6 ML capacity anti-spill jar to prevent leakage
- Color: Clear or green
- Diameter: 5mm or 6mm
- Kink-resistant tubing
- Latex-free
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Vernebler Maske	masque nébuliseur	máscara nebulizador	ネブライザーマスク


Simple Oxygen Masks



- Ref no:1.14.3.1
- Made from non-toxic PVC, clear and soft.
- Adjustable nose clip assures comfortable fit
- Adult and Pediatric sized available
- Color: transparent or green
- Diameter: 5mm or 6mm
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Kontakt xy Gen oder sks	Masques à oxygène	Las máscaras de oxígeno	酸素マスク

Oxygen Supply Tube



- Ref no:1.14.4.1
- Clear or green color
- Diameter: 5mm or 6mm
- Resist kinks and occlusion
- Length can be made as your request
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Kontakt xy Gen Tsube	Tube d'oxygène	tubo de oxígeno	酸素チューブ

Non-rebreathing Oxygen Masks



- Ref no:1.14.5.1
- Color: Clear or green
- Diameter: 5mm or 6mm
- Safety vent allows for entrainment of room air
- Adjustable nose clip assures comfortable fit
- 500ml-1000ml reservoir bag
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Kontakt xy Gen oder sks	Masques à oxygène	Las máscaras de oxígeno	酸素マスク

Venturi Mask



- Ref no: 1.14.6.1
- Safe, simple delivery of variable oxygen concentrations
- Features 6 color-coded diluters: 24%, 28%, 31%, 35%, 40% and 50%
- Includes adaptor for high-humidity entrainment
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Kontakt xy Gen oder sks	Masques à oxygène	Las máscaras de oxígeno	酸素マスク

Tracheostomy Mask



- Ref no. 1.14.7
- For tracheostomy and laryngectomy aerosol therapy. Tubing connector swivels 360 degrees and accepts 22mm I.D. corrugated tubing.
- Soft, clear vinyl mask features a tubing connector for tracheostomy and laryngectomy aerosol therapy.
- Adjustable elastic head strap ensures a comfortable fit.
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
Tracheotomie Mask	trachéotomie Mask	Máscara de traqueotomía	気管切開マスク

Humidifier bottle



- Ref no: 1.14.8
- Disposable, durable plastic with preset audible alarm at 4 -6psi pressure relief
- 250-500 mL usable volume
- Defined maximum and minimum water level lines
- Recessed nipple reduces possibility of breaking off outlet connection
- Packing: 1pc/polybag, 50pcs/ctn or as your required

			
anfeuchten Flasche	Humidifier bouteille	botella del humidificador	ウェットボトル

Oxygen bag



- Ref no: 1.14.9
- Main specifications : 30L 42L 50L
- That selected by the non-toxic, no harmful chemicals role of the single rubber-coated nylon material from a system of oxygen bags, main hospital for emergency transportation and family beds oxygen
- With clean, high-strength, easy to carry, safe, durable characteristics

			
Oxygen bag	Sac d'oxygène	bolsa de oxígeno	酸素バッグ

Aerosol Spacer



- Ref no: 1.14.10.
- Optimizes the delivery of MDI asthma medication.
- Compatible with most MDI (metered dose inhaler) actuators.
- Helps target medicine to your lungs.
- Helps keep you from spraying drug on the back of your throat.
- Anti Static plastic, Latex free
- With different sizes of mask for pediatric, child, adult use
- Packing: 1pc/box, 50pcs/ctn

			
Aerosol Spacer	aérosol Spacer	aerosol espaciador	エアゾールスプレーサー

3 ball breath exciser



- Ref no: 1.14.11
- Wide flow rate range of 600cc to 1200cc per second.
- Single patient use to reduce the risk of cross-contamination.
- Increase your inhaled lung volume and improve removing mucus or sputum.
- Different color-coded balls in three separate cylinders indicate the approximate flow rate passing through the device.

			
Atem exerciser	souffle exerciseur	ejercitador de respiración	息エクサ

PVC Manual Resuscitator



- Ref no:1.14.12
- All the best features at a price that fits your budget.
- Soft cushion mask allows for an effective seal on the patients face.
- The patient port swivels 360 degrees.
- Infant, pediatric and adult sizes available, with masks
- Reservoir style is available in bag or tubing.

			
Manuelle Beatmungsbeutel	insufflateur manuel	reanimación manual	マニュアル蘇生

Silicon Manual Resuscitator



- Ref no:1.14.13
- All the best features at a price that fits your budget.
- Soft cushion mask allows for an effective seal on the patients face.
- The patient port swivels 360 degrees.
- Infant, pediatric and adult sizes available, with masks
- Reservoir style is available in bag or tubing.

			
Manuelle Beatmungsbeutel	insufflateur manuel	reanimación manual	マニュアル蘇生