

Heat sealing flat reel



- Ref no. : 1.19.1.1
- Material: 55gsm PE film+60gsm paper
- Size: 50mm*200m ,100mm*200m can be customized
- Packing: 5-10rolls/ctn, can be customized
- Color option: blue

			
Heißsiegel flache Spule	La chaleur d'étanchéité bobine plate	carrete de sellado plana de calor	フラットリールヒートシール

Heat sealing gusseted reel



- Ref no. : 1.19.1.2
- Material: 55gsm PE film+60gsm paper
- Size: 50mm*200m
- Packing: 4rolls/ctn, can be customized
- Color option: blue or green

			
Heißsiegel gusseted reel	Bobine d'étanchéité thermique gousset	carrete de sellado por calor de fuelle	ヒートシールマチ付きリール

Self-sealing flat pouch



- Ref no. : 1.19.2
- Material: 55gsm PE film+60gsm paper
- Size: 70mm x230mm, 90mm x230mm, 135mm x255mm can be customized
- Packing: 200pcs/inner box,30 innerboxes/ctn
- Color option: blue



			
Selbstdicht flachen Beutel	Auto-étanchéité poche plate	Auto-sellado de la bolsa plana	自己シールフラットポーチ