

Scoop stretcher



- Ref.no.:4.3.1
- The stretcher can be divided into left and right halves by the clutches in both ends
- Adjustable length according to patient's height, easy to lock&unlock, and foldable
- It is made of high-strength aluminum alloy
- Can be used together with head immobilize
- Mainly used for ambulance, hospitals, emergency centers
- Accessories: 3pcs belts packed in a small bag
- Unfolded size(L*W*H) :214*42*7cm,Folded size:167*42*7cm,Load bearing:160kgs

			
Schaufeltrage	civière	camilla cuchara	スクープストレッチャー

Scoop stretcher



- Ref.no.4.3.2
- It has extra low profile for sliding the stretcher easily beneath the patient to minimize the risk of further injury
- Adjustable length according to patient's height, easy to lock&unlock, and foldable
- It is made of high-strength aluminum alloy
- Can be used together with head immobilize
- Accessories: 3pcs belts packed in a small bag
- Unfolded size(L*W*H) :214*42*7cm,Folded size:167*42*7cm,Load bearing:160kgs

			
Schaufeltrage	civière	camilla cuchara	スクープストレッチャー

Folding Stretcher



- Ref. no.:4.3.3
- It is made of high-strength aluminum alloy and oxford leather materials
- Optional color: blue, orange
- Optional fabric: oxford leather, orange PVC
- Unfolded size(L*W*H):185*52*4cm Folded size:93*52*8cm,Load bearing:160kgs

			
Faltstreckvorrichtung	Folding Stretcher	Camilla plegable	折りたたみストレッチャー

Folding Stretcher



- Ref.no.:4.3.4
- It is made of high-strength aluminum alloy and oxford leather materials
- Optional color: orange ,blue, camouflage, army green
- Optional fabric: oxford leather, orange PVC
- Unfolded size(L*W*H):221*53*11cm Folded
- Size:110*18*11cm,Load bearing:160kgs

			
Faltstreckvorrichtung	Folding Stretcher	Camilla plegable	折りたたみストレッチャー

Folding Stretcher



- Ref.no.:4.3.5
- It is made of high-strength aluminum alloy and oxford leather materials
- Optional color: orange ,blue, camouflage, army green
- Optional fabric: oxford leather, orange PVC
- Unfolded size(L*W*H):221*53*15cm Folded
- size:221*7*15cm,Load bearing:160kgs

			
Faltstreckvorrichtung	Folding Stretcher	Camilla plegable	折りたたみストレッチャー

Spine boards



- Ref.no.:4.3.6
- Material:PE•Runners keep the board raised for easy access lift
- Larger handholds designed to accommodate gloved hands
- Molded-in speed clip pins maximize strapping options
- Lightweight one-piece molded board
- 100% X-ray translucent feature gives users an extra clear view of the spinal area
- Compatible with most head immobilizers
- Easy to clean
- Packing:1pc/ctn

			
S Pine Bohrung rds	S Pine portait rds	tableros espinales	スパインボード

Patient mover



- Ref.no.4.3.7
- This compact mover is constructed from heavy-duty, non-woven materials and is able to carry up to 1,500 lbs.
- The fabric is impervious to fluids and the 14 handles allow for extra handlers when carrying larger patients.
- Packing:1pc/ctn

			
Pachien t Verschieben r	Pachien t Déplacer r	motor paciente	患者ムーバー

Extrication Device



- Ref.no.4.3.8
- Straps are color coded for ease of application.
- Can be used on pregnant and pediatric patients.
- Complete package includes head support, reusable head straps, and carrying bag.

			
Bergungsgerät	extrication Dispositif	dispositivo extricación	救出デバイス

Stair chair



- Ref.no.4.3.9
- The Stair Chair comes with three quick release buckle patient restraints,
- Easy to clean soft vinyl back seat and telescoping front handles
- Material:leather,aluminium,plastic
- Packing:1pc/ctn

			
Treppenstuhl	chaise d'escalier	S quieren r r chai	階段チェア