

Two Layer Folio Drape (Biplex Drape)




Ref no.: 1.9.1.1

- Material: PE film+ viscose or PP+PE or customized
- Size: 50*50cm, 45*75cm or customized.
- Packing: customized
- Style: disposable, sterile
- Color option: blue

			
zweischichtige foliant Drapierung	deux couches drapé folio	Para Aie r D Foio la Paix	2層フォリオドレーブ

Field Drape




Ref no.: 1.9.1.2

- Material: 16gsm paper+10gsm PE+16gsm paper.
- With or without hole.
- Use hot-pressing technology.
- Size: 45*65cm or customized
- Packing: customized
- Style: disposable ,sterile

			
Feld drapieren	champ drapé	campo drapeado	フィールドドレーブ

Paper+PE



Ref no.: 1.9.1.3

- Material: Paper+PE or 2 layer Paper+PE
- With or without hole (circle, Rhombus, Rectangle)
- With or without adhesive tape
- Color: white or blue
- Style: disposable ,sterile

			
Papierbeschichtung	Paper Coating	Recubrimiento de papel	紙加工

Three Layer Drape Sheet (Triplex Drape)




Ref no.: 1.9.1.4

- Material: Hydrophilic PP+ meltblown fabric +PE film or customized
- Size: 35*50cm, 50*50cm or customized
- Packing: customized
- Style: disposable, sterile
- Color option: blue

			
dreischichtige Drapierung Blatt	Trois couches Drapé Sheet	Th Rey Aie r D la Paix Shee t	三層ドレープシート

Arthroscopy Drapes (Knee Drape)




Ref no.: 1.9.1.5

- Material: 45gsm SMS with liquid collection pouch
- FDA Approved
- Size : 200*300cm or customized
- Packing: 20pcs/ctn
- Style: disposable, sterile
- color option: blue

			
Arthroskopie Drapierung (Knie Drapierung)	Trou de serviettes de l'articulation du genou	Toalla de la articulación de la rodilla	ひざ関節洞

Ophthalmic drape



Ref no.: 1.9.1.6

- Material: 40gsm SMS
- By CE/ISO Approved
- Packing: 20pcs/ctn
- Style: disposable, sterile
- Color option: blue

			
Augenheilkunde . Drapierung	Trou de tissu ophtalmique	Toalla de oftalmología	眼科洞

Split Drape (U Drape)



Ref no.:1.9.1.7

- Material:28gsm viscose+ 27gsm PE film
- Style: disposable, sterile
- high-protective
- Size : 80*100CM
- Packing: 1PC/POUCH, 100PCS/CTN

			
U - Drapierung	Trou de serviette en forme de U	Agujero en forma de U.	U 型ホール

Extremity Drape



Ref no.:1.9.1.8

- Material:45-50gsm SMS, SMMS
- With circular hole retractable film around incise area
- Size : 80*100CM
- Packing: 1PC/POUCH , 100PCS/CTN

			
Extremität Drapierung	La serviette de trou	Toalla de manos y pies	両手足の穴

Laparotomy Drape



Ref no.: 1.9.1.9

- Material:45gsm SMMS with fenestration incise film
- With adhesive tape around the fenestration
- Size : 80*100cm
- Packing: 1 set/pouch , 6 sets/ctn

			
Laparotomie Drapierung.	Trou de serviettes abdominales	Toalla abdominal	腹部穴

C-section Drape



Ref no.: 1.9.1.10

- Material: 45gsm SMMS, blue, with fluid collecting pouch, incise film.
- With triangular/round fenestrated incise area and circular fluid collection pouch and filtration sponge.
- Size: 260*300cm
- Packing: 1 pcs/pouch, 20 pouchs/ctn

			
C-Abschnitt Drapierung	La césarienne trou de serviette	La cesárea toalla	帝王切開洞

Angiography Drape



Ref no.: 1.9.1.11

- Material: glued with surgery film, water-absorbent material or outside membrane, no woven fabric,
- with Two holes or three holes.
- Size : customized
- Packing: 1pcs/pouch, 100pcs/ carton or customized

			
Angiografie Drapierung	Angiographie de trou de serviette	Angiografía toalla	血管造影

Cardiovascular Drape



Ref no.: 1.9.1.12

- Material: SMS medical angiography drape.
- With full chest and femoral incise area, adjustable tube holder and instrument pouches
- Size : customized
- Packing: 1pcs/pouch, 100pcs/ carton or customized

			
Kardiovaskulär Drapierung.	Trou de tissu cardio - vasculaire	Toalla cardiovascular	心血管心臟心臟洞

Half/Large and Utility Sheets



Ref no.: 1.9.1.13

- Dapes sheets are effective add-ons to other drapes to help create a sterile field and complete your procedural draping needs.
- Size : customized
- Packing: 1pcs/pouch, 100pcs/ carton or customized



üblich Blätter



Trou de serviette universel



En general la toalla



ユニバーサルホール

Disposable mayo stand Cover



Ref no.: 1.9.1.14

- Material: PE film+PP nonwoven
- Size : 78cmx145cm,80cmx140cm or Customized
- Packing: customized



Einmal Mayo Stand
Abdeckung



Mayo unique



Mayonesa simple



梅奥羽書

Table cover



Ref no.: 1.9.1.15

- Material: PE or SMS with or without reinforcement
- high quality
- Size : 150*190CM or customized
- Packing: 30PCS/CTN
- Non sterile



Wallpapers



Fonds d'écran



fondos de pantalla



壁紙

Disposable surgical urology drape




Ref no.: 1.9.1.16

- Material: 45-55gsm SMS or SMMS
- high quality
- Size : 150*240CM or customized
- Packing: 30PCS/CTN

			
Einmal chirurgisch Urologie Drapierung	Serviette de trou d'urologie	La Sección de la toalla de la cueva de urología	泌尿科洞

Utility drape




Ref no.: 1.9.1.17

- Material: Impervious PE +absorbent tissue paper
Or Impervious PE +Non woven material
- EO/Radiation
- Machine folding with fan
- Color: blue or green.
- Size : customized
- Packing: customized

			
üblich Drapierung.	General de deux couches de trou de serviette	El agujero de la capa de dos toallas	G M二段穴

Fenestrated Drape




Ref no.: 1.9.1.18

- Material: Impervious PE+absorbent tissue paper,
with Hole: oval, roundness or rectangle.
- Impervious PE+Non woven
- EO/Radiation
- Size : customized
- Packing: customized

			
Fenestrated Drapierung	Fenestrated Drapé	fenestrado cubrir	開穴巾

Utility drape with adhesive and hole



Ref no.: 1.9.1.19

- Material: Impervious PE +Viscose, hole: 8*8cm.
or Impervious PE+absorbent PP Diameter: 10cm
- with Medical double-sided tape
- Sterile
- Size : 50*50cm or customized
- Packing: customized

			
üblich Drapierung mit Klebstoff und fenester	General avec de la colle de la serviette.	Abrir agujeros con pegamento universal.	汎用ゴム開穴巾